REF NUS-001

Instructions for Use
Instrucciones de uso
Gebrauchsanweisung
Istruzioni per l'uso
Instruções de uso

HyperBranch Medical Technology, Inc.
800-12 Capitola Drive
Durham, NC 27713 USA

CE 0050

MedPass International Limited
Windsor House
Barnwood
Gloucester GL4 3RT
United Kingdom

REF NUS-001

Read Instructions for Use Prior to Using this Product.

Instructions for Use

Description:

The Adherus[®] Dural Sealant device packaging consists of two pouches with sterile contents and two sterile vial adapters. One pouch includes the precursor components for the preparation of a synthetic, absorbable hydrogel sealant and is labeled Adherus Dural Sealant. The other pouch includes the contents for the Spray Applicator for mixing and delivery of the precursor components to the target site and is labeled Spray Applicator Kit.

The pouch labeled Adherus Dural Sealant contains a foil pouch with two vials which contain the crosslinking components and a Tyvek pouch containing two 5 mL syringes with reconstitution solutions specific for each of the crosslinking components. Once the crosslinking components have been reconstituted with their respective reconstitution solutions, a polyethylene glycol (PEG) ester solution and a polyethyleneimine (PEI) solution (referred to as the 'green' and 'clear' precursors, respectively) are produced for subsequent assembly into the Spray Applicator.

The other pouch contains the Spray Applicator materials specific to mixing the precursors as they exit the tip of the Spray Applicator. When the solutions within the syringes are mixed together within the tip of the Spray Applicator, the precursors quickly crosslink to form the hydrogel sealant. The resulting hydrogel sealant is absorbed over approximately 90 days, sufficient time to allow for healing.

Indication:

Adherus Dural Sealant is intended for use as an adjunct to standard methods of dural repair, such as sutures, during neurosurgical intervention to provide watertight closure.

Contraindications:

- There are no known contraindications for this device.

Warnings:

- Adherus Dural Sealant is intended for use on the dura mater. The effectiveness of Adherus Dural Sealant on other tissue types has not been studied.
- The safety and effectiveness of Adherus Dural Sealant has not been studied in:
 - Patients with a known allergy to FD&C Blue #1 and/or FD&C Yellow #5.
 - Patients who are pregnant or lactating
 - Patients with severely altered renal or hepatic function.
- Do not use Adherus Dural Sealant if an active infection is present at the surgical site.
- Do not use Adherus Dural Sealant as a void filler in the spine, as post-operative hydrogel swelling may impinge on surrounding tissues.
- Adherus Dural Sealant is intended for use during neurosurgical procedures as an adjunct to standard methods of dural repair and should not be used as a replacement of standard methods of repair.

Cautions:

- The Adherus Dural Sealant pouches and blister packs are packaged sterile. Do not use if the pouches or blister packs have been damaged or opened. Do not re-sterilize.

- The Adherus Dural Sealant is intended for single patient use only. Discard opened and unused product.
- Do not use if the PEG powder is not free flowing.
- Use Adherus Dural Sealant within 2 hours of preparation of the green precursor.
- Prior to application of the Adherus hydrogel, ensure that adequate hemostasis has been achieved.
- Do not exceed approximately 2 mm of thickness. It may not be necessary to use the entire device.
- Avoid incidental application of Adherus hydrogel to tissue planes that will be subsequently approximated, such as muscle and skin.

Sterilization Method:

The Adherus Dural Sealant pouch contents and the vial adapters are sterilized with radiation. The Spray Applicator pouch contents are sterilized with ethylene oxide.

Storage:

Adherus Dural Sealant should be stored at or below 86 °F (30 °C).

How Supplied:

The components of the Adherus Dural Sealant package are:

- Tyvek Pouch Containing Sterile Adherus Polymer Kit (1)
 - ◆ Tyvek Pouch Containing Sterile Syringes (1)
 - Component A-2 (Green) PEG buffer in green plunger syringe with clear cap (1)
 - Component B-2 (Clear) PEI buffer in clear plunger syringe with clear cap (1)
 - ◆ Foil Pouch Containing Sterile Vials (1)
 - Component A-1 PEG ester powder in green capped vial (1)
 - Component B-1 PEI in silver capped vial (1)
- Tyvek Pouch Containing Sterile NUS001-2 Spray Applicator Kit (1)
 - ◆ Tyvek Pouch Containing Sterile (1)
 - Blending Connector (1)
 - Spray Tip (2)
 - ◆ Tyvek Pouch Containing Sterile (1)
 - Syringe Holder (1)
 - Plunger Cap (1)
- Vial Adapters (2) in Blister Pouches

DIRECTIONS FOR USE

The application procedure consists of three steps:

- 1) Preparing the Adherus Dural Sealant,
- 2) Assembling the Micromedics Spray Applicator system, and
- 3) Applying the Adherus Dural Sealant

STEP 1: PREPARING THE ADHERUS DURAL SEALANT

SIDE A Reconstitution

- 1) Using aseptic technique, introduce the contents of the large Tyvek peel pouch (foil peel pouch and a Tyvek peel pouch), labeled Adherus Dural Sealant, into the sterile field.
- 2) Open the foil pouch labeled Component A-1 and Component B-1 and open the Tyvek pouch labeled Component A-2 (Green) and Component B-2 (Clear) (see Figure 1.1). Set aside Components B-1 and B-2.

Figure 1: Adherus Dural Sealant pouch contents in the sterile field

- 3) Remove and discard the cap on the green syringe containing Component A-2 (syringe with a green plunger).
- 4) Outside the sterile field, remove the cover from the vial adapter package. Do not remove the vial adapter from the blister package.
- 5) Attach the sterile green syringe containing Component A-2 to the sterile vial adapter through the complementary luer-lock connections; use the blister pack to handle the adapter outside the sterile field (see Figure 2.1).
- 6) Remove and discard the blister package.
- 7) Seat the vial adapter with attached green syringe on the green vial containing Component A-1 by grasping the vial adapter and pushing down until the spike penetrates the septa and the adapter snaps in place (see Figure 2.2). Care should be taken not to depress the green syringe plunger during this process.
- 8) Once the vial adapter is seated, transfer Component A-2 to the green vial by slowly depressing the green plunger all the way in (see Figure 2.3). Retract the green syringe plunger to relieve the pressure within the glass vial (see Figure 2.4).
- 9) While holding the vial of Component A-1, shake and/or swirl the green solution vigorously for **at least one minute**. During this reconstitution step, air bubbles are introduced into the solution so that the solution will appear opaque even though the powder has dissolved.
- 10) Depress the green syringe plunger to remove all air within the syringe, invert the vial/syringe assembly, and withdraw the reconstituted Component A-1 from the vial into the green syringe (see Figure 2.5).
- 11) Remove the green syringe from the adapter.
- 12) Allow the air to dissipate from the green solution while preparing SIDE B. The solution should become transparent green with time.


Figure 2: SIDE-A Reconstitution

SIDE B Reconstitution

- 1) The glass vial containing Component B-1 appears empty but does in fact contain a small amount of liquid.
- 2) Remove and discard the cap on the syringe containing Component B-2.
- 3) Outside the sterile field, remove the cover from the vial adapter package. Do not remove the vial adapter from the blister package.
- 4) Attach the sterile syringe containing Component B-2 to the sterile vial adapter through the complementary luer-lock connections; use the blister pack to handle the adapter outside the sterile field (see Figure 3.1).
- 5) Remove and discard the blister package.
- 6) Seat the vial adapter with attached syringe on the vial containing Component B-1 by grasping the vial adapter and pushing down until the spike penetrates the septa and the adapter snaps in place (see Figure 3.2). Care should be taken not to depress the syringe plunger during this process.
- 7) Once the vial adapter is seated, transfer Component B-2 to the vial by slowly depressing the plunger all the way in (see Figure 3.3). Retract the syringe plunger to relieve the pressure within the glass vial (see Figure 3.4).
- 8) While holding the vial of Component B-1, shake and/or swirl the solution vigorously for **at least one minute**.
- 9) Depress the syringe plunger to remove all air within the syringe, invert the vial/syringe assembly, and withdraw the reconstituted Component B-1 from the vial into the syringe (see Figure 3.5).
- 10) Remove the syringe from the adapter.


Figure 3: SIDE-B Reconstitution

STEP 2: ASSEMBLING THE MICROMEDICS SPRAY APPLICATOR SYSTEM

- 1) Using aseptic technique, introduce the contents of the Spray Applicator Kit (two Tyvek pouches) into the sterile field (Figure 4.1).
- 2) Open the Tyvek pouch containing the Blending Connector and two Spray Tips and open the Tyvek pouch containing the syringe holder and the plunger cap (Figure 4.2).


Figure 4: Applicator Components

- 3) With the luer-lock tip for each syringe containing the Adherus SIDE-A and SIDE-B pointing up, tap the syringes to displace any air and then express the air out of the syringes.
- 4) Express the necessary amount of solution from the syringe containing a larger volume of solution so that the volume of solution within each syringe containing the Adherus SIDE-A and SIDE-B are equal (see Figure 5.1).


Figure 5: Applicator Assembly

- 5) Individually attach each syringe to the Micromedics Blending Connector via luer-lock connections (see Figure 5.2 and 5.3). Attach the syringe containing Adherus SIDE-B (clear solution) to the luer-lock connection designated on the applicator with an “0” (see Figure 5.3).
- 6) Each syringe, once in the threads, should be rotated clock-wise approximately 360° until the syringe flanges run perpendicular to the plane of the applicator (Syringe flanges are shown perpendicular to the applicator in Figure 5.4).
- 7) Attach the syringe holder onto the syringe barrels. The syringe flanges should easily fit into the groove in the middle of the syringe holder, located between the finger holds (see Figure 5.5).
- 8) Carefully attach the plunger cap to the plungers of both syringes without dispensing precursors through the blending connector (see Figure 5.6).
- 9) Assemble one of the included white Spray Tips to the Blending Connector by twisting firmly clockwise onto the luer-lock connection (see Figure 6.1). Two tips are provided.

STEP 3: APPLYING THE ADHERUS DURAL SEALANT

- 1) To ensure optimal performance, fluid outflow (CSF, blood, air) from the target site should be suspended during Adherus Dural Sealant application.
- 2) Hold the applicator approximately 2 to 4 cm from the target site. Apply firm even pressure to the center of the plunger cap to dispense Adherus. Applying a small amount of Adherus to a piece of gauze just off the target site and quickly moving to the target site during an application is recommended (see Figure 6.2).
- 3) Apply Adherus forcibly to enhance mixing, moving quickly along the target site. Care must be taken to ensure that both syringes are depressed simultaneously.
- 4) Continue applying the hydrogel until a thin coating (**approximately 1 – 2 mm**) is formed (see Figure 6.3).

Note: Gauging Thickness: Ensure that all suture knots are completely covered with hydrogel sealant. Typically, size 4-0, size 3-0 and size 2-0 sutures are used for dural closure. The smallest of these is size 4-0 which has a diameter of 0.15 mm to 0.2mm. A knot of size 4-0 suture will have at least four widths of suture or approximately 0.6 to 0.8 mm of thickness. Complete knot coverage ensures that the minimum thickness of application is achieved.

- 5) To stop delivery, release pressure on the syringe plungers. Note: The tip will plug.
- 6) If delivery of Adherus is interrupted and the Spray Tip is plugged, remove the plugged Spray Tip, wipe any hydrogel from the exposed tube ends, attach the second Spray Tip to the Blending Connector, and continue delivery.
- 7) Excess Adherus Dural Sealant beyond the edges of the dural margin may be removed with a Penfield probe, scissors or mechanical disruption.


Figure 6: Final Applicator Assembly and Application of Adherus Dural Sealant

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