Adherus[®]

AUTOSPRAY

REF NUS-006

Instructions for Use	4
Instrucciones de uso	14
Gebrauchsanweisung	25
Istruzioni per l'uso	36
Instruções de uso	47
Kullanma Talimatı	58

HyperBranch Medical Technology, Inc.
800-12 Capitola Drive
Durham, NC 27713 USA

CE 0050

MedPass International Limited
Windsor House
Barnwood
Gloucester GL4 3RT
United Kingdom

Explanation of Symbols / Explicación de los símbolos / Erläuterung der Symbole / Definizi dei simboli / Explicação dos símbolos / Sembollerin Açıklaması

The following table shows symbols used in this manual, on system components and on packaging labels. The table gives the corresponding meaning of each symbol. Pay attention to these symbols as they may relate to important safety precautions and warnings.

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Table of Contents

Device Description	5
Indication	6
Contraindications	6
Warnings	6
Cautions	6
Adherus AutoSpray Sterilization Method and Device Disposal	6
Environmental Operating Conditions	7
How Supplied	7
Pre-operative Procedures	7
Initial Inspection and Set Up	7
Adherus AutoSpray Operations	7
Reconstituting the Crosslinking Components	7
Turning the Adherus AutoSpray Device On	8
Adherus AutoSpray Non Replaceable Parts	8
Device and Treatment Site Preparation	9
Treatment Delivery	9
Stopping a Spray Application	9
Finishing a Delivery	10
Adherus AutoSpray Disposal	10
Troubleshooting Adherus AutoSpray	10
Contact Information for Further Assistance.....	10
Adherus AutoSpray Electrical Specifications	10
Guidance & Manufacturer’s Declaration – Electromagnetic Emissions	11
Guidance and Manufacturer’s Declaration – Electromagnetic Immunity for all Equipment and Systems	11
Guidance and Manufacturer’s Declaration – Electromagnetic Immunity	12
Recommended Separation Distances Between Portable and Mobile RF Communications and Adherus AutoSpray	13
General Technical Precautions	13

Caution: All users must read and understand the Adherus AutoSpray Instructions for Use, including indications, contraindications, operating instructions, warnings and precautions before performing any procedure. Failure to do so may result in injury to the patient or operator or cause damage to the Adherus AutoSpray device.

Device Description


Adherus AutoSpray is a sterile, single-use, electromechanical, battery operated, device with internal system components that provide air flow to aid in the delivery of a synthetic, absorbable, two-component hydrogel sealant system and allow delivery to be interrupted without clogging.

The device is supplied as a pre-assembled applicator and two separate glass vials, one of which is packaged within a foil pouch. The two glass vials contain either an activated polyethylene glycol (PEG) ester powder or a polyethyleneimine (PEI) dissolved in sterile water. The crosslinking components are reconstituted prior to use by their respective reconstitution buffers which are housed within the applicator. The resulting solutions mix within the applicator and quickly crosslink to form the hydrogel sealant soon after exiting the applicator tip.

Just before application, the Adherus AutoSpray applicator first must be turned on and primed. Then, while aiming at the treatment site, the syringe pusher assembly may be depressed to mix the two component system and deliver the resulting solution with the assistance of air flow. Once an application begins, it may be interrupted at any time by manually relieving the pressure on the syringe pusher assembly.

The delivered solution immediately crosslinks to form a hydrogel sealant which is absorbed over approximately 90 days, sufficient time to allow for healing.

The Adherus AutoSpray device is comprised of the following primary components:

Spray Nozzle

The spray nozzle thoroughly mixes the two sealant solutions and delivers the mixed solution to the target site through a tight spray pattern. The spray nozzle is attached to the system and is not removable.

Vial Adapter Hubs

The vial adapter hubs accept the vials containing the crosslinking components. During the reconstitution phase, the vials are attached to the vial adapter hubs. After reconstitution is completed, the hubs are removed from the Adherus AutoSpray device by rotating them counterclockwise. Removing the hubs opens the pathways for the solutions to flow through the nozzle.

On/Off Push Button Switch

The On/Off switch turns the battery operated air pump on and off. The device is shipped with the switch in the off position which isolates the air pump from the battery power source. Before removing the vial adapter hubs, push the ON/OFF switch to ON and the air pump will be engaged.

Battery Removal Door

A door on the underside of the housing which allows the operating room (OR) staff to remove the batteries for appropriate disposal at the end of the use if necessary. The battery door is glued shut. Use a flat instrument to pry the battery door open and remove the batteries for disposal.

Syringe Pusher Assembly

The syringe pusher assembly mechanically locks the two syringe plungers such that advancement of both syringe plungers occurs simultaneously.

Indication

Adherus AutoSpray is intended for use as an adjunct to standard methods of dural repair, such as sutures, during neurosurgical intervention to provide watertight closure.

Contraindications

- There are no known contraindications for this device.

Warnings

- Power is supplied by two AAA batteries.
- Adherus AutoSpray is intended for use on the dura mater. The effectiveness of Adherus AutoSpray on other tissue types has not been studied.
- The safety and effectiveness of Adherus AutoSpray has not been studied in:
 - Patients with a known allergy to FD&C Blue #1 and/or FD&C Yellow #5.
 - Patients who are pregnant or lactating
 - Patients with severely altered renal or hepatic function.
- Do not use Adherus AutoSpray if an active infection is present at the surgical site.

Cautions

- Adherus AutoSpray is packaged sterile. Do not use if the packaging or seals have been damaged or opened. Do not re-sterilize.
- Adherus AutoSpray is intended for single patient use only. Discard opened and unused product. Reuse of the device may result in cross contamination which may lead to injury, illness or death of the patient. Reuse may also exceed the 2 hour reconstitution time and be prevented by polymerization of the sealant within the applicator tip.
- Do not use if the PEG powder is not free flowing.
- Use Adherus AutoSpray within 2 hours of reconstituting the crosslinking components. When tested 8 hours after reconstitution, the resulting sealant swelled a small but statistically significant amount more when compared to samples reconstituted for 1, 2, or 4 hours.
- Prior to application of the Adherus hydrogel, ensure that fluid (cerebrospinal fluid, blood etc.) outflow has been suspended.
- Do not exceed approximately 2 mm of thickness. It may not be necessary to use the entire device.
- Avoid incidental application of the Adherus hydrogel to tissue planes that will be subsequently approximated, such as muscle and skin.

Adherus AutoSpray Sterilization Method and Device Disposal

The following table details the method of sterilization and disposal of the Adherus AutoSpray device. Because disposable components are provided sterile for initial use, do not use if the component packaging is open or damaged in any way.

Description	Part Number	Sterilization Method	Frequency of Use	Method of Disposal
Adherus AutoSpray	 NUS-006		Single Use Only 	Dispose of the Adherus AutoSpray device in a hospital approved biohazard container. Do not place in the trash! 

Environmental Operating Conditions

Storage Temperature	Store below 30°C
Operating Temperature	5°C to 35°C
Noise Level	≤60 db

How Supplied

The components of the Adherus AutoSpray package are:

- Adherus AutoSpray kit tray (1)
 - ◆ Pre-Assembled applicator (1)
 - ◆ Foil pouch containing glass vial (1)
 - Activated PEG powder in green capped vial (1)
 - ◆ PEI solution in silver capped vial (1)
 - ◆ Oxygen absorber (Do not remove the oxygen absorber from the tray)

Pre-operative Procedures

The following section provides instructions for receipt of the Adherus AutoSpray, initial inspection and set up. The user is advised to set up the Adherus AutoSpray as recommended in order to achieve safe and optimal use.

Initial Inspection and Set Up

Carefully inspect labels on the Adherus AutoSpray devices to verify the expiration date. If the system package is opened or damaged in any way, do not use the system and contact HyperBranch Medical Technology, Inc.

Caution: Keep all Adherus AutoSpray devices in shelf boxes and pouches until ready to use to protect them from damage or contamination.

Caution: Do not use the Adherus AutoSpray if beyond its expiration date.

Caution: Do not use the Adherus AutoSpray if the packaging is open, broken or compromised in any manner.

Caution: Do not use the Adherus AutoSpray in the presence of flammable anesthetics or flammable anesthetics with oxidants.

Caution: Keep the Adherus AutoSpray device away from strong magnetic fields to avoid possible interference of RF communication. Magnetic resonance equipment may generate interference with pump operations.

Caution: Medical electrical equipment requires special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the charts at the end of this users guide.

Caution: The Adherus AutoSpray device has a degree of protection against electric shock of applied part classified as a Type BF. This device is classified as internally powered medical electrical equipment.

Adherus AutoSpray Operations

The following sections give instructions for preparing and powering Adherus AutoSpray.

Reconstituting the Crosslinking Components

- 1) Using aseptic technique, introduce the contents of the large foil tear pouch labeled Adherus AutoSpray into the sterile field.
- 2) Remove the lid from the Adherus AutoSpray kit tray.
- 3) Remove the pre-assembled Adherus AutoSpray applicator and the adjacent vial containing the clear liquid from the tray.
- 4) With the device oriented such that the spray nozzle is pointing up, remove and discard the caps covering the spikes on the blue vial adapter hubs on the applicator (Fig. A-1).

Note: It may be necessary to use hemostats or a similar tool to remove the caps. A twisting motion while removing the cap is also recommended.

- 5) Seat vial containing PEI solution (**silver** capped vial) into the blue vial adapter hub of the applicator corresponding to the **white** syringe plunger by grasping the applicator and fully depressing the vial so that the spike penetrates the septa and the silver aluminum seal bottoms out in the vial adapter hub (Fig. A-2).
- 6) Remove the foil pouch containing the PEG ester powder from the tray.
- 7) Remove vial containing the PEG ester powder (**green** capped vial) from the foil pouch and seat vial into the blue vial adapter hub of the applicator corresponding to the **green** syringe plunger by grasping the applicator and fully depressing the vial so that the spike penetrates the septa and the green aluminum seal bottoms out in the vial adapter hub (Fig. A-2).
- 8) With the vials oriented upward, push the syringe pusher assembly to advance the two syringe plungers and transfer the liquids within the syringes into the vials. Continue to depress the syringe pusher assembly (Fig. A-3).

Note: Use two hands to depress the syringe pusher assembly if necessary.

- 9) With the syringe pusher assembly still depressed, gently shake the device intermittently until the powder is completely dissolved (Fig. A-3).

Note: The PEG ester powder will continue to dissolve in between periods of gentle shaking; overly aggressive shaking can cause excessive entrapped air bubbles.

- 10) Allow the syringe pusher assembly to rebound backwards and push forward a second time to ensure both solutions are homogeneous.


- 11) Allow the syringe pusher assembly to rebound backwards again and pull the syringe pusher assembly backwards until all of the reconstituted ingredients are evacuated from the vials (overhead view in Fig. A-4).

Note: If the On/Off Push Button is accidentally engaged during the reconstitution phase, the device may be turned off with no detrimental effects.

Caution: Use must be completed within two hours of reconstitution for optimal results.

Turning the Adherus AutoSpray Device On

- 12) Push the ON/Off button on the bottom of the device to “activate” the system prior to hydrogel spray application (Fig. B-1).

Note: The Adherus AutoSpray applicator should be activated prior to removal of the vial adapter hubs, just prior to spray application, to avoid inadvertently clogging the nozzle.

Adherus AutoSpray Non Replaceable Parts

Adherus AutoSpray is powered by two internally sealed AAA batteries that are non-replaceable.

Caution: Adherus AutoSpray is a single use sterile product. Do not re-sterilize.

Caution: Portable and mobile RF communications equipment can affect medical electrical equipment.

Device and Treatment Site Preparation

13) If not already on, push the ON/Off button to turn on the Adherus AutoSpray applicator (Fig. B-1).

! **Note:** The internal air pump is automatically activated by pushing the ON/Off button, providing an audible sound that the system has been activated. If the audible sound from the pump is not present, discontinue use of the current Adherus AutoSpray unit and follow instructions in the Troubleshooting section.

14) Prior to treatment delivery, orient the device so that the vials are pointed upwards and remove and discard both vial adapter hubs by rotating the blue hubs counterclockwise until the hubs are ejected from the device (overhead view in Fig. B-2).

15) For optimal tissue adherence at the treatment site, ensure that 2 to 3 millimeter margins around the defect are clear of blood, hemostatic agents or other loose tissues and that cerebrospinal fluid outflow is minimized.

Treatment Delivery

16) To apply the prepared product, if not already on, push the ON/Off button in to turn on the device.

17) Orient the device such that the spray nozzle is pointed up and hold a piece of gauze approximately 5 to 10 centimeters above the nozzle (Fig. B-3).


18) Apply firm even pressure to the center of the syringe pusher assembly until product begins to spray out of the nozzle (Fig. B-3).

19) Once a green sealant begins to form on the piece of gauze, stop depressing the syringe pusher assembly. Adherus AutoSpray is now ready for the treatment delivery.

20) While aiming at the treatment site and holding the device nozzle approximately 2 to 4 centimeters away, apply firm even pressure to the center of the syringe pusher assembly to dispense the mixed solution (Fig. B-4).

21) Continue applying the Adherus sealant system until a thin coating (**approximately 1 – 2 mm**) is formed.

! **Note:** Gauging Thickness: Ensure that all suture knots are completely covered with hydrogel sealant. Typically, size 4-0, size 3-0 and size 2-0 sutures are used for dural closure. The smallest of these is size 4-0 which has a diameter of 0.15 mm to 0.2mm. A knot of size 4-0 suture will have at least four widths of suture or approximately 0.6 to 0.8 mm of thickness. Complete knot coverage ensures that the minimum thickness of application is achieved.

Stopping a Spray Application

22) Once a spray application begins, the application may be interrupted at any time during a treatment by reducing pressure on the pusher assembly. However, do not turn off the pump until the final sealant application has been completed.

! **Note:** If excess hydrogel forms on the spray nozzle during an application, it may be removed by gently wiping with a piece of gauze.

Finishing a Delivery

- 23) When the final spray application is complete, turn off the device, open the battery door to break the internal seal and remove the batteries for appropriate disposal.
- 24) Excess Adherus Dural Sealant beyond the edges of the dural margin may be removed with a Penfield probe, scissors or mechanical disruption.

Adherus AutoSpray Disposal

The following table relays important information regarding disposal of the Adherus AutoSpray components.

Device	Disposal Instructions
Adherus AutoSpray	Adherus AutoSpray consists of a plastic enclosure and electrical components. This device must be taken to separate collection at the product end of life. Do not dispose of this product as unsorted municipal waste. 

Troubleshooting Adherus AutoSpray

The following potential failure modes for Adherus AutoSpray have been identified. Each of the failure modes includes a built-in safeguard (as listed below) to control for the side effects of such malfunctions.

Potential Failure Modes	Safeguard
1. Air-pump not activated	1. On/Off switch not depressed
2. Spray nozzle blocked	2. Do not remove vial adapter hubs prior to turning on air pump, wipe off nozzle after use.

 **Caution:** If a problem occurs with Adherus AutoSpray, identify the symptom then attempt to resolve the problem as indicated in the Instructions for Use. If the problem cannot be resolved, under NO condition should application continue in the hope that the problem will be remedied by itself.

To avoid complications, follow exactly the Adherus AutoSpray Troubleshooting Flow Chart in the above table. During a procedure, the user's efforts to troubleshoot Adherus AutoSpray should be strictly limited to the steps described above. Discontinue use of the current Adherus AutoSpray unit **immediately** when indicated, open a new Adherus AutoSpray unit, and contact HyperBranch Medical Technology, Inc. following the procedure.

Contact Information for Further Assistance

In the event that the device becomes inoperable and cannot be remedied using the troubleshooting guidelines above, contact:

HyperBranch Medical Technology, Inc.
800-12 Capitola Drive
Durham, North Carolina 27713-4410
919-433-3325

Have the system lot number (located on the Adherus AutoSpray labeling) available for reference when contacting HyperBranch Medical Technology, Inc.

Adherus AutoSpray Electrical Specifications

The Adherus AutoSpray meets the following standards:

- IEC 60601-1 3rd edition Electrical safety standard requirements
- IEC 60601-1-2 3rd edition Immunity requirements

Guidance & Manufacturer's Declaration – Electromagnetic Emissions

Adherus AutoSpray Sealant is intended for use in the electromagnetic environment specified below. The customer or user of Adherus AutoSpray should assure that it is only used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	Adherus AutoSpray uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Adherus AutoSpray is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations /Flicker emissions IEC 61000-3-3	Not applicable	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for all Equipment and Systems

Adherus AutoSpray is a battery operated system and is intended for use in the electromagnetic environment specified below. The customer or user of Adherus AutoSpray should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	Not Applicable Battery Operated System with no Digital Circuitry	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not Applicable Battery Operated System	Main power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	Not Applicable Battery Operated System	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) for 5 sec	Not Applicable Battery Operated System	Main power quality should be that of a typical commercial or hospital environment. If the user of the equipment requires continued operation during power main interruptions, it is recommended that the equipment be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not Applicable Battery Operated System with no Magnetically Sensitive Components	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

*NOTE: U_T is the a.c. main voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

Adherus AutoSpray is intended for use in the electromagnetic environment specified below. The customer or user of Adherus AutoSpray should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80MHz to 2.5GHz</p>	<p>Not Applicable Battery Operated System with no Digital Circuitry</p> <p>Not Applicable Battery Operated System with no Digital Circuitry</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of Adherus AutoSpray, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended Separation Distance</p> <p>$d = 1.2\sqrt{P}$</p> <p>$d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.5GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicated theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Adherus AutoSpray is used exceeds the applicable RF compliance level above, the Adherus AutoSpray device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Adherus AutoSpray device.</p> <p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less the 3 V/m</p>			

Recommended Separation Distances Between Portable and Mobile RF Communications and Adherus AutoSpray

Adherus AutoSpray is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Adherus AutoSpray device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Adherus AutoSpray device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Caution: The use of accessories, transducers and cables other than those specified by the Manufacturer may result in increased Emissions or decreased Immunity of the Adherus AutoSpray device.

Caution: Adherus AutoSpray should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the Adherus AutoSpray device should be observed to verify normal operation in the configuration in which it will be used.

General Technical Precautions

Caution: Do not remove system covers. The manufacturer or an authorized service representative only must perform installation, assembly, service and maintenance. Adherus AutoSpray is a one-time use disposable product.

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